



## Welcome to the BITalino Family

Here's some useful information to get you started with:

[www.bitalino.com](http://www.bitalino.com) is the place to go to:



- download [Software](#) and [APIs](#) for data acquisition
- access [Documentation](#) to get things going
- have the latest info about [BITalino](#)




## Experience your BITalino in 4 easy steps.\*

- 1 Turn on BITalino
- 2 Pair BITalino via Bluetooth setup (PIN: 1234)
- 3 Add the MAC address in the OpenSignals software settings
- 4 Record your signals!

\* FAQ and detailed setup instructions can be found in the user manual, which can be download from the Help Screen inside the OpenSignals software or from [bitalino.com](http://bitalino.com).

## Here are a few examples of what you can do...



Wearable Drone Controller



Robotic Prosthetic Hand



Holter-like ECG device



Excitement-reaction Room

## Have any questions?

Visit [forum.bitalino.com](http://forum.bitalino.com) and share your **thoughts** and **projects**.

We'd also love to get your **feedback**: